
ES-720D

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7 Push button "P" (7). In order to visualise the preset temperature, press and release it. It is possible to change the value using the UP and DOWN buttons. At the end of the programming phase the read in is automatically done. If the bottom is pushed for a too long time it come to seek on programming: in this case, not change any parameter, the item come to seek normal after few seconds. / Pousser le bouton "P". Afin de visualizer la température pré réglée, presser et relâchez-le. Il est possible de changer la valeur en utilisant les bouton UP et DOWN. A la fin de la phase de programmation la lecture est effectuée automatiquement. Si le bouton est poussé pour un temps trop long il va chercher sur la programmation: dans ce cas, ne pas changer les paramètres, l'élément va chercher la normale après quelques secondes.

DIGITAL THERMOSTATIC UNIT

The Digital Thermostatic Unit **ES- 720 D** is equipped with an helical mechanical stirrer and heating element of 800 W. Electronic thermoregulator to regulate the room temperature at 80°C, resolution ± 1 °C visualised on a 2 figures digital display unit.

Protection device for insufficient fluid level inside the tank. Outer casing in powdered painted stainless steel Aisi 304, acids resistant.

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