



ESBD-M7DP Disposable Video Laryngoscope with 7 Inch Monitor



1

Single use blades, reducing the risk of cross-infection, reprocessing costs, time and storage.

2

Anti-Fog Design

3

Separate 7inch LCD Monitor, to keep doctor away from patients secretions.

4

More Available Space with 12mm

5

Blade work with both ESBD-M3 portable monitor and ESBD-M7 7inch Monitor



Model Number: **ESBD-M7DP**

| Order number | ESBD-M7 | ESDP-MILO | ESDP-MAC1 | ESDP-MAC2 | ESDP-MAC3 | ESDP-MAC4 |
|------------------|-----------------|-----------|-----------|-------------|-----------|-------------|
| Description | 7in LCD Monitor | Mil0 | MAC1 | MAC2 | MAC3 | MAC4 |
| Application | LCD Monitor | Neonate | Pediatric | Small Adult | Adult | Large Adult |
| Blade Length(cm) | / | 8.31 | 10.15 | 12.11 | 16.08 | 17.5 |

Technical Data

| Parts | Description | Technical Index |
|-----------------------------|---------------------------------------|--|
| Monitor | Display | 3 inch LCD |
| | Resolution Ratio | 640*480(RGB) |
| | Viewing angle | >60° |
| | Power | <2W |
| | Display front and back rotation angle | 0°~130° |
| | Monitor Thickness | 12mm |
| | Weight | <90g |
| Camera | Illuminance | >150 LUX |
| | Resolution Ratio | 1280*720 px |
| Image/Video function | Image/ Video function | Yes |
| | Output | USB output, Easy to storage and establish image /video |
| | Memory type | Micro SD card 8GB |
| | Video/Photo Play back | Yes |
| Battery | Battery Type | Rechargeable lithium battery |
| | Capacity | 730mAh |
| | Cycle life of battery | >300 times |
| | Battery working time | >75min |
| | Charging time | <90min |
| Blades | Weight | < 80g |



ESSE3 srl, Via Garibaldi 30
 14022 Castelnuovo D.B. (AT)
 Tel +39 011 99 27 706
 Fax +39 011 99 27 506
 e-mail: danielec@esse3-medical.com
 web: www.esse3-medical.com

