



Endoscope Cabinet



The inner tendon of the storage cabinet is made of stainless steel materials. The soft mirror cabinet is designed with a endoscopic suspension system made of transparent acrylic. Prevent collisions with each other. And the lower parts are lifted to meet the demand for endoscopic storage of different sizes.

Model	Product	Type	Specification(mm)	Features
FBS-200	Respiratory	single door	170x55x50	UV lamp, Biopsy hook,tangential ventilation turbines
FBS-400		double door	170x90x50	
GC-200	Gastro	single door	180x55x50	
GC-400		double door	190x90x50	
CC-200	Colono	single door	200x55x45	
CC-400		double door	210x90x50	



ESSE3 srl, Via Garibaldi 30
 14022 Castelnuovo D.B. (AT)
 Tel +39 011 99 27 706
 Fax +39 011 99 27 506
 e-mail: danielec@esse3-medical.com
 web: www.esse3.dreamgest.com

