



## Shaver System



### Code: 367.004 Shaver System

- Speed range 500-3500 RPM.
- Variable modes forward, reverse & oscillating modes.
- Foot switch. easy foot control in variable blade modes (forward, reverse & oscillating).
- Shaver hand piece with card fully autoclavable.
- Suction control inbuilt in hand piece.
- Power supply 230V (AC) / 50 Hz



## Pulse Lavage System



**Code: 367.005**

### **Pulse Lavage System**

Irrigation pressure electronically controlled.

- Autoclavable and reuseable tubes of medical grade silicon.
- Complete sterility is maintained.
- Easy loading of tubes in the machine.
- Efficient peristaltic pump.
- Foot control operated machine.
- Pre - set features.
- Four types of irrigaton possible.
- Inbuilt fuse and voltage control rectifier.
- Powersupply 230v AC 50 Hz  $\pm$  10 %
- Plastic Body.