



ES-KE180

Electric Dialysis Chair



BLOOD BANK

ES-KE180

Electric Dialysis Chair

Technical Configuration

- * Hand Controller-----1pc
- * Power Line-----1pc
- * 3" Silent Caster-----1set
- * Small dining table-----1pc
- * With CPR
- * Four 3" castors with break

Main Technical Parameters

Seat size: 1900mm × 580mm (L*W)
 Seat height adjustable: 550mm-750mm
 Back size: 770*580mm(L*W)
 Seat size: 470*580mm(L*W)
 Leg size: 380*540mm(L*W)
 Foot board size: 410*200mm(L*W)
 Seat height: 550-750mm
 Mattress thickness: 6cm
 Back rest lifting: 0 ° -70 °
 Leg rest lifting: 20°-75°
 Armrest swing out: 0 ° -270 °
 Armrest size: 500 × 160mm(L*W)
 Armrest height adjustable: 100mm

