




7 < GYfJYg

Hand pieces with
NSK connection



DENTISTRY



7 < 'GYf]Yg

Hand Pieces



With NSK connection

Cod. CH-MUR



- | | |
|------------------|----------------------------------------------------------------------|
| • DESCRIPTION | Standard head high-speed handpiece, Wrench type, Single spray system |
| • AIR PRESSURE | 180-250 KPa (1.8-2.5 kgf / cm ²) |
| • WATER PRESSURE | 198Kpa (2 kg) |
| • ROTATION SPEED | 370,000- 420,000 rpm |
| • HEAD HEIGHT | 13.14 mm |
| • HEAD DIAMETRE | Ø 11.2 mm |
| • CHUCK TYPE | Wrench |
| • BUR APPLICABLE | Ø1.59mm ~ Ø1.6mm x 21mm~23mm (diameter x length) |
| • NOISE | ≤70dB |

Cod. CH-MIC




- | | |
|------------------|---------------------------------------------------------------------|
| • DESCRIPTION | E type low speed handpiece, 1:1 Ratio, with external cooling system |
| • AIR PRESSURE | 245~392Kpa (2.5~4.0 kgf/cm ²) |
| • WATER PRESSURE | 198Kpa (2 kg) |
| • ROTATION SPEED | 22,000~27,000 rpm |
| • BUR APPLICABLE | Ø2.334~2.355 mm |
| • NOISE | ≤70dB |

Cod. CH-CAB



- | | |
|------------------|---------------------------------------------------------------------|
| • DESCRIPTION | E type low speed handpiece, 1:1 Ratio, with external cooling system |
| • AIR PRESSURE | 245~392Kpa (2.5~4.0 kgf/cm ²) |
| • WATER PRESSURE | 198Kpa (2 kg) |
| • ROTATION SPEED | 22,000~27,000 rpm |
| • BUR APPLICABLE | Ø 2.334~2.355 mm |
| • NOISE | ≤70dB |

Cod. CH-CAV



- | | |
|------------------|----------------------------------------------------------------------------------------------------------|
| • DESCRIPTION | 4:1 reduction angle handpiece, with external cooling system connect to E Type motor, precise gear drive. |
| • CHUCK TYPE | Latch |
| • BUR APPLICABLE | Ø2.35 mm |
| • NOISE | ≤70dB |



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

