



# ESE-TROLLEY

Trolley for ESE line  
ultrasound systems



IMAGING SYSTEMS



# ES9!HF C@@9M

'Hfc`YmZcf`9G9`]bY`i`hfUgci bX`gnghYa g



### Aluminium Alloy Material

- Light weight
- Robust & Resistant
- Less cost for transportation
- Long life circle to meet the complex emergency environment

### Tilted Table Design

- Better view angle
- Comfortable to operate the system on the trolley

### Gas Lifted Table

- Easy to handle
- Transformed according to your requirements whenever you sit or stand for operation

### Simple Assembling

It's easy for you to attach or detach the portable ultrasound system with the trolley. It's flexibility makes you move fast when you are busy with your work



4 probe holders



Extended 2 or 3 probe connectors (optional)



Extended USB ports



Physical keyboard bracket (optional)



Extended VGA/HDMI video port (optional)



**Battery pack**  
It can be the backup of the battery inside the system-lasts for 5 hours working time



Video printer bracket



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

