



# SONOTRAX

Ultrasonic Pocket Doppler



IMAGING SYSTEMS



- Easy- to- use
- Compact design
- Interchangeable waterproof probe
- Built-in audio recorder
- Large backlight LCD display
- Accurate FHR detection with high-fidelity sound effect
- Support waterproof probes for both obstetrical and vascular usage
- Auto, Average and Manual working modes
- Various power supply solutions
- Automatic power-off to save power
- Support 2MHz,3MHz,4MHz,5MHz,8MHz waterproof probe

**Display:** 45x25mm LCD display

**Heart rate range:** 50-210bpm

**Resolution:** 1bpm

**Audio sampling frequency:** 4KHz / 8KHz

**Recording length:** 480s / 240s

**Audio output power:** 0.5W

**Auto shut-off:** 1 minute

**Size:** 32x85x138mm D\*W\*H

**Weight:** 290gr including one battery

### SonoTrax Lite&Vascular

Interchangeable waterproof probes, alkaline battery

### SonoTrax Basic

Interchangeable waterproof probes, LCD display, alkaline battery

### SonoTrax Basic A

Interchangeable waterproof probes, LCD display, backlight, three working modes, alkaline battery

### SonoTrax Pro

Interchangeable waterproof probes, LCD display, backlight, three working modes, built-in audio recorder, alkaline battery

### SonoTrax II

Interchangeable waterproof probes, LCD display, backlight, three working modes, Li-ion rechargeable battery

### SonoTrax II Pro

Interchangeable waterproof probes, LCD display, backlight, three working modes, built-in audio recorder, Li-ion rechargeable battery



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

