



# ESFD-680

Fetal Doppler



INTENSIVE CARES



# ESFD-680

## Features

1. Durable and economic with volume wheel
2. Color TFT display or LCD display optional
3. Two interchangeable probes, High sensitive, anti-shock, waterproof, anti-fall
4. Built-in speaker, Headphone audio output
5. Automatic shut down after 1 minute
6. Battery status indicator and low power indication

## Specifications

Ultrasonic frequency: 2MHz+3% and 3MHz +3%  
 Ultrasonic intensity: <10mW/cm<sup>2</sup>

Power supply: 2 Common AA batteries  
 Continuous working time: > 10 hours

Working temperature: 10°C ~ 40°C  
 Dimension: 120mm\*110mm\*40mm

### Packaging:

Normal packaging: 100 Units/CTN  
 Carton size: 71.5x35x59.5cm

GW: 26.3KG

### Unit packaging:

Box size: 16.5x13.8x5.8cm

GW: 0.25KG (Not including battery)



ESSE3 srl, Via Garibaldi 30  
 14022 Castelnuovo D.B. (AT)  
 Tel +39 011 99 27 706  
 Fax +39 011 99 27 506  
 e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
 web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

