



Cod. **ZB-60**

Zoom Range	0.8X / 6.5X
Zoom Ratio	1/8
Eyepiece	WF10x/22mm
Objective	Zoom
Viewing head	Binocular
Illumination	LED

- Binocular head inclined at 20° with interpupillary distance (55-75mm) and Diopter adjustment.
- Eyepiece WF10/22mm whit ring rubber.
- Focusing system whit coaxial coarse.
- Working Distance of 78mm.
- LED reflected/transmitted illumination with brightness adjustment.
- Weight 5Kg.
- 1 White base disk.
- User Manual.





Cod. ZB-60

Portée du zoom	0.8X / 6.5X
Taux de zoom	1/8
Oculaire	WF10x/22mm
Objectif	Zoom
Tête de visualisation	Binoculaire
Illumination	LED

- Tête binoculaire inclinée à 20° avec réglage de la distance interpupillaire (55-75mm) et réglage dioptrique.
- Oculaire WF10/22mm avec anneau en caoutchouc.
- Système de mise au point coaxiale grossière.
- Distance de travail de 78mm.
- Éclairage à lumière réfléchie/transmise LED avec réglage de la luminosité.
- Poids 5Kg.
- 1 Disque base blanc.
- Manuel d'utilisation.

CE