




ESBDF Series

Combined refrigerator
and freezer



LABORATORY



Features

Name: Combined refrigerator and freezer

Model: ESBDF-25V450E / ESBDF-25V450F

Outside cabinet: PCM color plate

Outside door: Epoxy polyester powder spray coated steel plate

Condenser: Wire tube condenser

Evaporator: Refrigerator: Tube and plate evaporator

Freezer: Wire tube evaporator

Insulation: Rigid polyurethane foam filling

Compressor: Fully enclosed

Thermostat: Microcomputer control

Temperature sensor: NTC

Alarm system: High/low temperature alarm, door ajar alarm, power failure alarm, low battery alarm, sensor error alarm.

Specifications

Model	Ambient temperature	Climate type	Refrigerant	Rated voltage (V ~)	Rated frequency (Hz)	Storage temperature (°C)	Net capacity (L)	Rated current (A)	Net Weight (Kg)	Exterior dimension (W*D*H) (mm)
ESBDF-25 V450E	16~32°C	N	R600a/30g R600a/65g	220	50	Up: 2~8 Down: -10~-26	450	2.11	144	810×735×1960
ESBDF-25 V450F	16~32°C	N	R600a/30g R290/67g	220	50	Up: 2~8 Down: -10~-40	450	3.05	146	810×735×1960



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

