



ES-RM80 Series

Lab Tube Roller Mixer



LABORATORY



Description

ES-RM80 provides a gentle but highly efficient rocking and rolling motion, ideal for mixing blood samples, viscous substances and liquid-solid suspensions.

Used in variety applications, prevention of blood coagulation and immune precipitation etc. and not destroy blood cells.

Features

- Digital speed control and LED display
- Continuous or timer operation
- Excellent performance, adjustable speed and quiet running
- Tube in 8mm-200mm with a variety of applications
- The up/down movement ensures complete mixing of the sample
- Apply to 4°C cold chamber or 60°C incubator with relative humidity of 80%

Specification

Model	ES-RM80	ES-RM80D
Speed range	10rpm-80rpm	
Rocking motion	Rolling	Rolling/Swinging
Amplitude	/	25mm
Numbers of rollers	6	
Roller length	280mm	
Max. load capacity	4.2kg	
Timed range	1min~99h59m	
Voltage	AC 100~230V, 50/60Hz	
Power	30W	
Fuse	250V, 1A, Φ 5x20	
Dimensions (W × D × H)	(mm) 230x440x120mm	
Weight	3.8 Kg	



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

