



Pneumatic Tourniquets



LITTLE ACCESSORIES



Pneumatic Tourniquets

1 circuit or 2 circuits

Electro-mechanical pneumatic tourniquet - 1 circuit

Cod. G10705

Easy to use :

- 1 regulated pressure circuit
- Controlled by microprocessor, a guarantee of security and accuracy
- Setting pressure and timer
- Setting range : 0-600 mmHg
- Accuracy of ± 5 mmHg
- Timing and pressure alarms

Optional :

- Mobile stand on casters
- Arm and thigh cuffs - one bladder



Electro-mechanical pneumatic tourniquet - 2 circuits

Cod. G10706

Easy to use :

- 2 regulated pressure circuits
- Controlled by microprocessor, a guarantee of security and accuracy
- Setting pressure and timer
- Setting range : 0-600 mmHg
- Accuracy of ± 5 mmHg
- Timing and pressure alarms

Optional :

- Mobile stand on casters
- Arm and thigh cuffs - two bladders



Electronic pneumatic tourniquet - 1 circuit

Cod. G10803

Easy to use :

- Touch screen interface : 120mm x 90 mm
- 1 regulated pressure circuit
- Setting range : 0-600 mmHg
- Controlled by microprocessor
- Battery
- Compact

Traceability :

- Ease of maintenance
- All fuses are externally accessible
- Removable battery via an access hatch

Security :

- Pressure alarm and timer alarm (audible and visual)
- Battery alarm (audible and visual)

Optional :

- Printer
- Mobile trolley on casters
- Arm and thigh cuffs - 1 bladder



Electronic pneumatic tourniquet - 2 circuits

Cod. G10903

Easy to use :

- Touch screen interface : 120mm x 90 mm
- 2 regulated pressure circuits
- Setting range : 0-600 mmHg
- Controlled by microprocessor
- Battery
- Compact

Traceability :

- Ease of maintenance
- All fuses are externally accessible
- Removable battery via an access hatch

Security :

- Pressure alarm and timer alarm (audible and visual)
- Battery alarm (audible and visual)

Optional :

- Printer
- Mobile trolley on casters
- Arm and thigh cuffs - 2 bladders





Optional accessories



A10701

Mobile trolley on casters

for G10705, G10706, G10803 and G10903



G13200

Imprimante

pour G10803 et G10903

Brassards et Cuissards Optionnels



Arm and thigh cuffs - 1 bladder
for G10705 and G10803

Arm and thigh cuffs - 2 bladders
for G10706 and G10903

| | |
|----------------|-----------------------------------|
| GBDM101 | Straight cuff - Infant size - Arm |
| GBDM102 | Straight cuff - Child Size - Arm |
| GBDM103 | Straight cuff - Adult S - Arm |
| GBDM104 | Straight cuff - Adult M - Arm |
| GBCM105 | Straight cuff - Adult L - Arm |
| GCDM103 | Straight cuff - Adult S - Thigh |
| GCDM104 | Straight cuff - Adult M - Thigh |
| GCDM105 | Straight cuff - Adult L - Thigh |
| GCDM106 | Straight cuff - Adult XL - Thigh |
| GCDM107 | Straight cuff - Adult XXL - Thigh |
| GCCM104 | Conical cuff - Adult M - Arm |
| GCCM105 | Conical cuff - Adult L - Arm |
| GCCM106 | Conical cuff - Adult XL - Arm |
| GCCM107 | Conical cuff - Adult XXL - Arm |

| | |
|----------------|-----------------------------------|
| GBDM201 | Straight cuff - Infant size - Arm |
| GBDM202 | Straight cuff - Child Size - Arm |
| GBDM203 | Straight cuff - Adult S - Arm |
| GBDM204 | Straight cuff - Adult M - Arm |
| GBCM205 | Straight cuff - Adult L - Arm |
| GCDM203 | Straight cuff - Adult S - Thigh |
| GCDM204 | Straight cuff - Adult M - Thigh |
| GCDM205 | Straight cuff - Adult L - Thigh |
| GCDM206 | Straight cuff - Adult XL - Thigh |
| GCDM207 | Straight cuff - Adult XXL - Thigh |
| GCCM204 | Conical cuff - Adult M - Arm |
| GCCM205 | Conical cuff - Adult L - Arm |
| GCCM206 | Conical cuff - Adult XL - Arm |
| GCCM207 | Conical cuff - Adult XXL - Arm |



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

