



ESF2-359/00

Gynaecological couch
with variable height



MEDICAL FURNITURE



ESF2-539/00

Gynaecological couch with variable height



Available colors



Gynaecological couch with variable height

- Adjustable head and leg sections in any position.
- Upholstered top covered with washable white plastic material.
- Frame made in sheet steel epoxy painted with adjustable feet.
- Variable height from 56 to 93 cm.
- Removable armrests.
- Complete with S.S. steel bowl \varnothing 32 cm., and two adjustable knee crutches.

Cod. **ESF2-539/00** Standard version with electrical variable height

Cod. **ESF2-539/10** Version with hydraulic variable height

Optional wheels



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

