



ESF2-2862/00

Delivery bed



MEDICAL FURNITURE



ESF2-2862/00

Delivery bed



Cod. ESF2-2862/00



Optional



Optional

Technical specification:

- 3-Section foam mattress with high density covered with PU leather, hygienic and durable.
- Square tube epoxy coated steel frame with anti rust treatment.
- Leg section mounted on silent sliding rails under the seat.
- Adjustments of back rest by ratchet.
- 80 cm fixed height.
- Fixing Clamp for anatomical knee crutches.
- Large stainless steel sliding tray for liquid collection under the seat.
- Central part (Seat) with trendelenburg position for easy and comfortable delivering.
- Demountable sturdy construction.
- All section are separated and studied with safety distance and safety block to avoid injuries due to improper use or emergency situation.
- Safe working load 230 Kg.
- Trendelenburg 0 – 15°
- Length 180 cm
- Width 80 cm
- Height 80 cm



Included accessories:

- Mattress Set
- Knee Crutches
- Padded shoulder rests
- Hand grips
- Stainless steel tray

Optional accessories:

- Wheels
- Side rails



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

