



ES-CACF Line

Dirty Linen Trolleys



MEDICAL FURNITURE



ES-CACF Line

Dirty Linen Trolleys



TECHNICAL SPECIFICATIONS

Cod.	ES-CACF203CC	ES-CACF203CC110
Overall dimensions	1090x690x1515h mm	1410x730x1515h mm
Internal dimensions	1000x600x1270h mm	1320x640x1220h mm
Volume	720 L	1100 L
Weight	44 kg	51 kg
Wheels	Diam. 160 mm Diam. 200 mm	
Opening	Equipped with nr.4 wheels in grey anti-mark rubber diam.160 or 200 (2 fix and 2 swivel), disposition at the corners.	
Opening	Tip-downdoor and lid for a total closure	
Ribs & Push handle	Horizontal reinforcement ribs each 120 mm approx. Horizontal push handle fixed at one head.	
Material	Made in light alloy plate, anodised after trolley finished, Base chassis in light alloy profiles	

Optional accessories

Cod. PARA	Anti-mark PVC bumper in light alloy profile, along the whole lower trolley perimeter profile.
Cod. FREINS	Bumpers for 2 wheels



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

