



ESB!%5





Medical bed - 4 sections



MEDICAL FURNITURE




Cod. ES6-%



Control button for fold the lateral rails



Folded side rails



Opened crank



Bed end with 4 cranks

- The backrest, the legrest, the Trendelenburg / Anti-Trendelenburg positions and the height are adjustable by four hand cranks
- The platform of the bed is made of perforated steel plate with anti-corrosion treatment
- Epoxy painted steel frame with anti-rust treatment
- The ABS ends panels are easily removable
- 4 silent rotating casters, two with brakes
- Easily foldable security side rails in aluminium alloy
- 4 anti-shock bumpers on the four corners of the bed
- Telescopic I.V. pole on the side of the bed

Technical features

- Trendelenburg / Anti-Trendelenburg position : -12/+12°
- Backrest tilting angle : 1-80°
- Legrest tilting angle : 1-40°
- Adjustable height range : 250mm
- Casters diameter: 125 mm
- Max. patient weight : 200 Kg
- Sleeping surface dimensions : 190 x 55 x 25 cm
- Overall dimensions : 208 x 98 x 45~70 cm