



# STAINLESS STEEL CARTS

## CHARIOTS EN ACIER INOX

**Cod. 816026**



**Cod. 816027**



Cart entirely made from glazed stainless steel AISI 304. The cart is particularly robust, handy, with removable shelves with raised edges to avoid spilling.



*Chariots entièrement réalisés en acier inox AISI 304 satiné. Chariots particulièrement robustes, maniables avec des plateaux amovibles à rebord.*

## TECHNICAL SPECIFICATIONS CARACTERISTIQUES TECHNIQUES

Cod	SIZE (mm)	TOTAL HEIGHT (mm)	N° SHELVES	SWIVEL CASTERS Ø (mm)
	<i>TAILLE (mm)</i>	<i>HAUTEUR TOTALE (mm)</i>	<i>N° PLATEAUX</i>	<i>ROUES PIVOTANTES Ø (mm)</i>
<b>816026</b>	650 x 450	810	2	100
<b>816027</b>	650 x 450	810	3	100
<b>816028</b>	800 x 500	835	2	100
<b>816029</b>	800 x 500	835	3	100
<b>816030</b>	850 x 550	835	2	100
<b>816031</b>	850 x 550	835	3	100
<b>816032</b>	1070 x 600	870	2	125
<b>816033</b>	1070 x 600	870	3	125

### OPTIONAL ACCESSORIES ACCESSOIRES OPTIONNELS

	S.S. RAISED EDGE	S. S.BOTTLE HOLDER	WASTE BAG HOLDER	S.S. TRAY WITH SUPPORT	SET OF 4 BUMPERS	SET OF 4 ANTISTATIC CASTORS
	<i>PETITE BARRE</i>	<i>PORTE FLACONS</i>	<i>PORTE-SACS</i>	<i>CUVETTE EN INOX AVEC SUPPORT</i>	<i>SET 4 GALETSPARE-CHOCES SET 4</i>	<i>SET 4 ROUES ANTISTATIQUE S</i>
<b>Cod.</b>	<b>817025</b>	<b>817022</b>	<b>817023</b>	<b>817020</b>	<b>817024</b>	<b>817027</b>



**ESSE3** srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

