



# ESK-B041-3

Emergency & Recovery  
Trolley



MEDICAL FURNITURE



# ESK-B041-3

## Emergency & Recovery Trolley

### **DESCRIPTION:**

- Size(L\*W\*H): 2100\*750\*600-920mm
- Equiped with two USA hydraulic pumps
- Electrostatic spray steel frame and ABS base
- Foldable aluminium alloy guardrail
- With 100mm waterproof PU maffress
- 8 inch central controlled caster
- Pedal control elevation and
- trendelenburg/reverse trendelenburg located on both sides of the table
- X-ray translucent platform
- 5th wheel to control the direction of moving

### **FUNCTIONS:**

Back—rest lifting: 0°-82 ° (±5°)

Leg—rest lifting: 0°-30 ° (±5°)

Reversed trendelenburg: 0—18°

Trendelenburg: 0—18°

Height adjustment: 615-935mm

### **ACCESSORIES:**

- With oxygen cylinder holder and tray for placement of x—rays cassette
- With 1 piece SS telescoping I V stand and previsions
- With rubber bumpers on the four corners
- Loading bearing: 220kg



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

