



ESF1-800/xx

Baby Crib



ESF1-800/00 BABY CRIB

Trolley manufactured from chrome plated steel tube mounted on four casters Ø 60 mm.
Crib manufactured from Plexiglas, with aeration holes and with position TR/RTR.

ESF1-800/10 BABY CRIB

Trolley manufactured from stainless steel tube mounted on four casters Ø 60 mm.
Crib manufactured from Plexiglas, with aeration holes and with position TR/RTR.



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

