



ES-KM1202EL

C-arm compatible semi-electric
operating table



OPERATING THEATRE

ES-KM1202EL

C-arm compatible semi-electric operating table



Electric
movement
up/down

SPECIAL FEATURES:

- Lateral square column
- 304 Stainless Steel Fitting
- Detachable head/leg/pelvic section
- In built kidney bridge
- Five Sectional Radio – Translucent Top.
- Base & Cover made of attractive, easy to clean Stainless Steel.
- Sophisticated mechanics provide smooth, step less articulation of the table Top for precise patient positioning.
- Backup battery (OPTIONAL)



Trendelenburg



Reverse
Trendelenburg



Flex Position



Reflex Position



Lateral Tilt
Left & Right



Height Adjustment
Up & Down

TECHNICAL FEATURES

Length of Table	2000 mm
Width of Table Top	520 mm
Height (without Mattress)	Min. 750 to max. 1000 mm
Trendelenburg	+/- 25°
Lateral Tilt	+/- 20°
Back Section	+/- 70°
Head Section	+/- 60°
Leg Section	0°-90°
Floor Locking System	Manual

STANDARD ACCESSORIES



Arm Rest



Shoulder
Support



Leg Support



Anesthetic
Screen



Clamps



Mattress

STANDARD ACCESSORIES:

- Autoclavable mattress
- Pair of arm supports
- Anesthetic screen
- Pair of shoulder supports
- Pair of perineal (lateral) supports
- Pair of leg supports
- Set of clamps
- Wired remote control