



ESHE-205-3A

Ophthalmic electric
operating table



OPERATING THEATRE

ESHE-205-3A

Ophthalmic electric operating table



FEATURES:

This table is suitable for examination, diagnosing, surgical operation and ophthalmology under the microscope. It designs super-low position stable structure. The table has ophthalmology special accessories for EENT.

TECHNICAL FEATURES

Dimensions	Table top length	2030 mm
	Table top width	600mm
Functions	Height for table top (Electric)	580 -780 mm
	Electric tilt back and forth	+15°/ -35°
	Head board elevator adjustment (Electric)	0-120mm
Standard accessories	Anesthetic screen fram for ENT	1pc
	Arm supports for doctor	2pcs
	Arm supports for patient	2pcs
	Manual Controller	1pc
Optional accessories	Electric Operating Chair Electric elevator 550-750mm\Foot switch\it is movable and can be locked	
Materials	Stainless Steel base device, the mattress is made of sponge with high quality leather cover	