



YZ20P5

Operating Microscope



OPHTHALMOLOGY / ENT



Description

The product is a simple coaxial illumination operation microscope for one man.

It is portable, compact and applied in mobile hospital

Characteristics:

All the lenses are multi-coated, mildewproof and anti-reflective.

Optional objectives: F250/f300/f350f/f400.

Column stand can be reflitted to holding stand.

It can also be customized according to various operation requirements.

Specifications

Eyepiece Magnification: 12.5×

Objective Lens: f=200

Working Distance: 190mm

Magnifications for Main Microscope: 5.3×,8×,12×

Diameter of Field: F37mm,F25mm,F16.7mm

Diopter Adjustment: ±7D

Pupil Distance: 50mm~80mm

Maximal Resolution: 100LP/mm

Illumination Source: 12V/100W,medical cold reflection halogen lamp



ESSE3 Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

