



# ESR-300A

Ophthalmic unit



OPHTHALMOLOGY / ENT



# ESR-300A

## Ophtalmic unit

### Specification

#### **Mechanical Parameters**

Left/Right angle of swing arm:  $\pm 30^\circ$

Angle of rotation for swing arm right-angle rack:  $360^\circ$

Angle of rotation for Table:  $90^\circ$

Slide-board Slippage Range: 500mm

Lifting height for chair: 150~200mm

Loading for chair: 70~150kgs

#### **Electrical Parameters**

Input voltage: 110V/220V AC 50HZ~60HZ

Engine power: 55W 24VDC 40W 30VDC

Light luminaire: 18W 250V

Fuse: 5A 250V AC

Power without load: 1W

Pressure test: 2KV

Main wire: 10A 250V

Working mode: S1

Working temperature:  $-30^\circ \sim +50^\circ$

Relative humidity: less than 90%

Force of pressing key: 1~5N

Cooling mode: Nature cooling

#### **Installation Dimension**

Figure size: 1060mm×1000mm×950mm

Weight: 210kgs



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

