



Hand-held Slit Lamp Microscope



Specifications:

- Eyepiece: 10X, 16X
- Total magnification: 10X, 16X



- Diopter adjustment range : $\pm 7D$
- Work distance : 80mm
- Slit Width: 0~10mm, Continuously adjustable
- Filters: Green, cobalt blue, color temperature transformation angle
- Slit Angle: Horizontal circumference $\pm 30^\circ$
- Slit Illumination: White LED
- Input Voltage: AC 100~240V 60~50Hz
- Input power: 18VA
- Package dimensions: 195×105×230mm
- Gross weight: About 750g (including battery)



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com





Microscope lampe à fente Portable



Caractéristiques:

- Oculaires: 10X, 16X
- Grossissement total: 10X, 16X



- Champ de réglage de la dioptrie : $\pm 7D$
- Distance de travail : 80mm
- Largeur de la fente: 0 ~ 10 mm, réglable en continu
- Filtres: vert, bleu cobalt, couleur de la température de l'angle de transformation
- Angle de fente: $\pm 30^\circ$ horizontal circonférence
- Système d'illumination: lampe LED blanc
- Tension d'entrée: AC 100~240V 60~50Hz
- Puissance d'entrée: 18VA
- Dimensions de l'emballage: 195×105×230mm
- Poids brut: 750g (compris la batterie)



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

