

# AC-600 GASTRIC-THORACIC DRAINAGE PUMP

## FEATURES & FUNCTIONS :

1. Adopt drum closed twin pump design. No noise, no bad smell, no spark, safety and long life.
2. Digital timer adopts micro computer programming design and large display for easy reading and simple operation.
3. The output pressure of pump can be instantly known from the pressure test push button.
4. The new watch type pressure gauge which can clear show the reading replaces general water manometer.
5. The vacuum bottle is made of impact resistant material, so it can be sterilized in high temperature and high pressure.
6. Germ killing filter can prevent germs intruding from outside to secure safety.
7. One large article box at the rear of instrument is used for storing the leading tube, catheter, etc.
8. Stainless steel skeleton for easy maintenance, diagonal PU antistatic brake wheels for easy fixing.
9. Control panel has safe protective cover.

## STANDARD ACCESSORIES :

1. Silicon suction tube (120 cm) .....
2. Filter (big) .....
3. Fuse .....
4. Two-way adapter .....

## SPECIFICATIONS :

Voltage	AC 120V or 240V
Power consumption	20W or 10W
Negative pressure regulator	0~1500 mm / H <sub>2</sub> O
Max. volume	2000c.c. x 2 pcs
Safe device	Level float control
Machine size	L40 x W40 x H83 (cm)
Packing size	L46 x W46 x H89 (cm)
Net weight	17.8 KGS
CE NO.	DF-500



Via Garibaldi 30 – 14022 Castelnuovo Don Bosco (AT) Italy

TEL: +39 011 992 77 06 FAX: +39 011 992 75 06

HTTP: [www.4ci.com](http://www.4ci.com) E-MAIL: [info@4ci.com](mailto:info@4ci.com)

