

AG-600 OBSTETRIC DELIVERY & OPERATING TABLE

FEATURES:

1. MULTI FUNCTION AS OBSTETRIC DELIVERY & OPERATION TABLE
2. THE TABLETOP FRAME, CYLINDER COVER, BASE COVER, AND ACCESSORIES ARE MADE OF STAINLESS STEEL 304 WHICH IS DURABLE AND EASY TO BE CLEANED & MAINTAINED.
3. THE 60CM CLEARANCE OFF THE GROUND PROVIDES EASY AND SAFE ACCESS TO THE PREGNANT WOMEN. THE 30CM HOIST FURTHER MEETS THE REQUIREMENT OF THE DOCTOR IN TERM OF THE HEIGHT
4. THE SPECIAL DESIGN OF TABLE BASE OFFER THE PRACTITIONER WITH GREATER SPACE FOR ACCESSIBILITY AND COMFORT DURING OPERATION
5. THE 57CM WIDE HEAVY DUTY TABLE TOP OFFER A COMFORT LYING AND SITTING POSITION FOR THE PREGNANT WOMEN. IT CAN BE AN OPERATION TABLE IMMEDIATELY WITHOUT MOVING THE PATIENT IF NECESSARY
6. PU FOAM LEG CRUTCHES OFFER SAFE AND COMFORTABLE FOR THE USER
7. FULLY ELECTRIC HYDRAULIC FOOT CONTROL SYSTEM WHICH IS WATER PROOF AND WITH 24V CONTROL SWITCH



AG-600 OBSTETRIC DELIVERY & OPERATING TABLE



STANDARD ACCESSORIES:

1. PILLOW:	1 PC
2. ARM BOARDS:	2 PCS
3. MATTRESS:	1 PC
4. WATER PROOF SEAT COVER:	1 PC
5. WAIST SIDE HANDLES:	2 PCS
6. LEG CRUTCHES:	2 SETS
7. AUXILIARY TABLE:	1 PC
8. SLIDING PAN (WITH FILTER):	1 PC
9. FOOT CONTROLLER:	1 SET
10. POWER CORD:	1PC

SPECIFICATIONS:

MODEL	AG-600
TABLETOP	57 X 124CM
TABLETOP WITH AUXILIARY TABLE	57 X 185CM
ELEVATION	62 - 92 CM
TRENDELENBURG	+30 - -10
BACK UP/DOWN	+75 - -10
POWER SUPPLY	AC110V OR 220V, 50/60HZ

PS: AG-601 IS THE SAME AS AG-600 BUT WITH ELECTRIC CONTROL HYDRAULIC BRAKE.



ESSE 3 - Via Garibaldi 30 - 14022 Castelnuovo DB (AT) Italy

TEL: +39 011 99 27 506

FAX: +39 011 99 27 506

HTTP: WWW.4ci.com

E-MAIL: info@4ci.com

