

ES – EOVX 2

Vacuum Extractor electrical

“MOBINOX” SERIE

- CLASS I
- TYPE B
- USE : CONTINUOUS OR BY DOUBLE PNEUMATIC FOOT
- HIGH FLOW – HIGH VACUUM
- POWER SUPPLY 230 V – 50 Hz
- CONSUMPTION 240 W
- CONFORMITY EN 60601.1 - ISO 10079.1 - 93/42/CEE - CEI 64-4
- MAX VACUUM 85 kPa
- MAX LIMITED FLOW 50 l/m'
- MAINTENANCE FREE MOTOR , MAX FLOW 70 l/m' AT 0,85 BAR, PROVIDED OF SAFETY OVER-HEAT DEVICE
- 18/10 STAINLESS STEEL FRAME
- ANTISTATIC WHEELS AT BASE
- 2 WAY SELECTOR SWITCH
- 2 AUTOCLAVABLE POLYCARBONATE JARS, EACH 3,2 LITERS, PROVIDED OF LIDS WITH SAFETY TRAP , STOP VACUUM DEVICE AND ANTIBACTERIAL FILTERS
- VACUUM MANOMETER DIAMETER 63 MM, MAX – 1,0 BAR (-100 Kpa)
- TRIPLE MICROMETRIC ADJUSTABLE VACUUM SWITCH FOR CONTINUOUS , PNEUMATIC FOOT, OFF
- SHELF FOR ACCESSORIES
- STERILISABLE PROBE GLASS HOLDER
- SIZE cm. 53x36x90 h
- WEIGHT Kg 24



ESSE 3 Via Garibaldi 30
14002 Castelnuovo D.B.(AT) – I
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

