

**SURGICAL SUCTION EQUIPMENT FOR
LIPID SUCTION APPLICATION**

Model **ES - LIPO 1** (code 0300101)
"MOBINOX" Series



Practical
Functional

Ready for any suction operation

Sturdy

(No special maintenance required)

General features :

- 18/10 stainless steel **stand** with rounded corners
- Antistatic and antiblock **rotating wheels**
- Antistatic rubber **bumper**
- **Collection container support** with rubber base and adjustable stainless steel springs
- **3 l. Polycarbonate collection container** which can be sterilized. The cover is fitted with an overflow valve. (Containers of different types and capacities can be used).
- **Antibacterial filter** – 24 hour continuous operation
- **Expansion chamber** to protect the suction circuit, with overflow valve and quick coupling for disconnecting the group.
- **Vacuum gauge** Ø 63 mm, scale : -100 / 0 kPa; -76 / 0 cm Hg
- **Bipolar switch** : Continuous – Off – Foot switch
- **Micrometric needle regulator**
- Maintenance free **suction pump**, membrane type, flow rate 70 l./min., max. vacuum level 0,85 bar = 85 kPa. With thermal protection. Single phase motor.
- **Noiseless operation**
- **Pneumatic foot switch**, with start/stop functions
- **Convenient tray** for holding instruments
- **Pirex glass test-tube** for holding tubing
- **Silicone tube** Ø 10 x 16 or at your choice, size to be stated in the order.
- **Double funnel connection**
- **Dimensions** : cm. 40x36x90 h
- **Weight** : Kg 20

CLASS I



Type B

OPERATION MODE : *CONTINUOUS*

High Flow – High Vacuum

Power supply 230 ~ 50 Hz
VA 200

FUSE T 2 A Ø 5x20



IN CONFORMITY WITH EN 60601.1
ISO 10079.1
93/42/CEE

Options :

Polycarbonate test tubes

IN USE: VACUUM -80 kPa
FLOW 40 l/m'

³ The equipment can also be used for surgical application if equipped with appropriate accessories.



ESSE 3 Via Garibaldi 30
14002 Castelnuovo D.B.(AT) – I
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

