

## SURGICAL SUCTION EQUIPMENT FOR LIPID SUCTION APPLICATION

Mod. **ES - LIPO 2** (cod. 0300102 )  
"MOBINOX" Series



Practical  
Functional  
Ready for any suction operation  
Sturdy  
(No special maintenance required)

**General features :**

- 18/10 stainless steel **stand** with rounded corners
- Antistatic and antiblock **rotating wheels**
- Antistatic rubber **bumper**
- **Collection container supports** with rubber base and adjustable stainless steel springs.
- Collection container **selector**
- **Two 3 l. capacity Polycarbonate collection containers** which can be sterilized. The cover is fitted with an overflow valve. (Containers of different types and capacities can be used).
- **Two Antibacterial filters** – 24 hour continuous operation
- **Two Expansion chambers** to protect the suction circuit, with overflow valve and quick coupling for disconnecting the group.
- **Vacuum gauge** Ø 63 mm, scale : -100 / 0 kPa; -76 / 0 cm Hg
- **Bipolar switch** : Continuous – Off – Foot switch
- Collection container **manual selector**
- **Micrometric needle regulator**
- Maintenance free **suction pump**, membrane type, flow rate 70 l./min., max. vacuum level 0,85 bar = 85 kPa. With thermal protection. Single phase motor.
- **Noiseless operation**
- **Pneumatic double foot switch**, with start, stop and vacuum level programming functions
- **Convenient tray** for holding instruments
- **Pirex glass test-tubes** for holding tubing
- Stabilized silicone **connecting tube** Ø 10 x 16
- **Double funnel connection**
- **Dimensions** : cm. 58x36x87 h
- **Weight** : Kg 24

CLASS I



TYPE B  
OPERATION MODE : *CONTINUOUS*  
**High Flow – High Vacuum**

Power supply 230 ~ 50 Hz  
VA 200

FUSE T 2 A 5x20 



IN CONFORMITY WITH EN 60601.1  
ISO 10079.1  
93/42/CEE

<sup>4</sup> The equipment can also be used for surgical application if equipped with appropriate accessories.

IN USE: VACUUM -80 kPa  
FLOW 40 l/m'

**Options**

**Polycarbonate test tubes**



ESSE 3 Via Garibaldi 30  
14002 Castelnuovo D.B.(AT) – I  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.4ci.com](http://www.4ci.com)

