

ES - SURGICAL ASPIRATORS



SUCTION ES 1

- Electric suction unit on trolley, one bottle
- Frame in stainless steel
- Antistatic wheels
- Free maintenance motor pump
- Autoclavable MAKROLON bottle , cap. 2 lt
(available capacity 4 lt)
- Max vacuum 950 mbar
- Adjustable vacuum 0-950 mbar
- Max flow 80 lt/min
- Control manometer 0-1000 mbar
- Security double trap
- Probe holder
- CEI 62/5 et IEC 601 norms
- Dimensions 65 x 35 x 125 cm , Weight 37 kg

Equipped with :

- Aspiration tube and probe
- STOP-VAC handle (optional)
- Waterproof pneumatic foot (optional)

ES 1T : top table suction unit

Same technical features of above , but top table version

- Carrying bar
- Dimensions 65 x 35 x 95 cm , Weight 31 kg



SUCTION UNIT ES 2

- Electric suction unit on trolley, double bottle
- Frame in stainless steel
- Antistatic wheels
- Free maintenance motor
- Double autoclavable MAKROLON bottle, capacity 2+2 lt (available capacity 4 lt)
- Max vacuum 950 mbar
- Adjustable vacuum 0-950 mbar
- Max flow 80 lt/min
- Control manometer 0-1000 mbar
- Security double trap
- Probe holder
- CEI 62/5 et IEC 601 norms
- Dimensions 85 x 35 x 125 cm , Weight 43 kg

Equipped with :

- Aspiration tube and probe
- STOP-VAC handle (optional)
- Waterproof pneumatic foot (optional)

DISTRIBUTED BY :



ESSE 3 Via Garibaldi 30
14002 Castelnuovo D.B.(AT) - I
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com





ES 2T MOBILE SUCTION APPARATUS

- Electric suction unit on trolley, double bottle
- Frame in stainless steel
- Antistatic wheels
- Free maintenance motor pump
- Autoclavable MAKROLON bottle , cap. 2 lt
(available capacity 3 and 4 lt)
- Max vacuum 850 mbar
- Adjustable vacuum 0-850 mbar
- Max flow 60 lt/min
- Control manometer 0-900 mbar
- Security double trap
- Probe holder
- CEI 62/5 et IEC 601 norms
- Power supply 220 V (optional built-in rechargeable battery)
- Dimensions 40 x 35 x70 cm , Weight 10 kg

Supplied with all accessories



ES-3 MOBILE SURGERY-ANAESTHESIST SUCTION UNIT

- Electric suction unit on trolley, triple bottle
- Independent adjustment
 - 2 bottles for surgery capacity 1t 3
 - 1 bottle for anaesthetist capacity 1t 1
 - (available capacity 3 and 4 lt)
- Frame in stainless steel
- Antistatic wheels
- Free maintenance motor
- Max vacuum 950 mbar
- Adjustable vacuum 0-950 mbar
- Max flow 90 lt/min
- Control manometer 0-1000 mbar
- Security double trap
- Probe holder
- CEI 62/5 et IEC 601 norms
- Dimensions 85 x 35 x 125 cm , Weight 57 kg

Supplied

- Hoses and probes
- footswitch

DISTRIBUTED BY :



ESSE 3 Via Garibaldi 30
14002 Castelnuovo D.B.(AT) – I
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

