

ES-350VB VERTICAL STEAM STERILIZER



FULLY AUTOMATIC
CONTROL SYSTEM
WITH DRY VACUUM



ES-350VB VERTICAL STEAM STERILIZER

ESSE 3 ES-350 VB vertical type steam sterilizers with fully automatic control system, and dry (post) function, have space-saving, easy operation, durable, large capacity features.

Designed for clinics and laboratories with versatile sterilizing needs and where objects are to be sterilized such as: test tubes, glassware, liquids, surgical instruments, garments, human waste products etc.

FEATURES:

1. Automatic pressure control from heating, sterilization, exhausting, drying to finish, and included with wrapped/unwrapped and liquid sterilization mode.
2. With wide opening (350mm) and shorter of chamber that making easy to operating sterilizing subject.
3. Sterilization temperature is adjustable from 121°C to 134°C on temperature/pressure regulator.
4. Sterilization time and drying time is adjustable from 0~60minutes on timer.
5. The timer starts only after the pre-set temperature is reached to ensure the complete sterilization.
6. Illuminated pilot lamps display status and a buzzer announces the completion and cuts off the heating power source automatically.
7. The pressure gauge make ideal require reading and indicate high precision and consistent results.
8. Chamber and door is made of SUS 316, and door holder is made of stainless steel SUS 304 with polished surface for long lasting service and easy cleaning.
9. Four heavy-duty castors are mounted for convenient movement.

MODEL	ES-350VB
OVERALL	W65 x D51 x H 108cm
CHAMBER DIAMETER / DEPTH	350 mm / 600 mm
CHAMBER SIZE	LITERS
STERILIZATION TEMPERATURE	121C -134 C
LIQUID STERILIZATION MODE	YES
DRY FUNCTION	YES
POWER CONSUMPTION	3.0KW
VOLTAGE/FREQUENCY	220V/50-60HZ
WEIGHT	About 94 KGS
SAFETY DEVICES	
OVER TEMPERATURE PROTECTION	YES
ELECTRIC OVERLOAD PROTECTION	YES
OVER PRE	YES
EMERGENCY EXHAUST SWITCH	YES
SAFETY VALVE	YES
OVER HEAT/TEMP. INDICATION LIGHT	YES

STANDARD ACCESSORIES:

1. 1. BOTTOM SHELF x 1PC
2. STAINLESS STEEL WIRE BASKET 330cm x H27cm x 2PC
3. OPERATION MANUAL x 1PC



ESSE 3 Via Garibaldi 30
 14002 Castelnuovo D.B.(AT) – I
 tel +39 011 99 27 706
 fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web : www.4ci.com

