



EST - 700 SERIE

OPHTHALMIC & PLASTIC OPERATING TABLE



HIGH QUALITY

ECONOMIC

STABLE



EST - 700 SERIE

OPHTHALMIC & PLASTIC OPERATING TABLE

FEATURES:

1. EST-700 SERIE IS HIGH QUALITY , AND WITH REAL ECONOMIC BENEFIT TABLES. IT IS ESPECIALLY SUITABLE FOR OPHTHALMIC AND PLASTIC OPERATING TABLES.
2. EST-700 SERIE IS DESIGNED AND MADE IN HIGH RIGID STRUCTURE WITH TERRIFIC STABLE PERFORMANCE.
3. EST-700 SERIE IS ACTIVATED BY ELECTRIC-HYDRAULIC SYSTEM, WITH WATER PROTECTED FOOT CONTROLLER.
4. TABLE FRAME, BASE, COLUMN, COVER, AND OTHER IMPORTANT UNITS ARE MADE OF STAINLESS STEEL FOR RUST-PROOF AND EASY CLEANING.
5. THE MATTRESS IS MOVABLE, SO THE HAND REST FOR SURGEON IS VERY EASY TO INSERT AND TO SET UP.
6. FOR EST-702, THERE IS A BACK ANGLE ADJUSTMENT FUNCTION. IT CAN MAKE PROPER ADJUSTMENT FOR POSITION IN OPERATIONS. IT ALSO OFFER A GREAT ASSISTANCE TO SOME WEAK PEOPLE TO GET OFF AFTER THE OPERATION.



EST-701



EST-702

STANDARD ACCESSORIES:

1. FOOT CONTROLLER: 1 SE
2. POWER CORD: 1 PC

OPTIONAL ACCESSORIES:

- 1 HAND REST OF HEAD SECTION FOR DOCTOR



EST-704

SPECIFICATIONS:

MODEL	AT-701	AT-702	AT-704
TABLETOP	600 X 1900 MM	600 X 1900 MM	600 X 1900 MM
ELEVATION	550-850 MM	550-850 MM	550-850 MM
BACK SECTION	NA	0-70 °	NA
POWER SUPPLY	AC110V OR 220V, 50/60HZ	AC110V OR 220V, 50/60HZ	AC110V OR 220V, 50/60HZ



ESSE 3 Via Garibaldi 30
 14002 Castelnuovo D.B.(AT) - I
 tel +39 011 99 27 706
 fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web : www.4ci.com

