

ESSE 3 ENDOSCOPE

PURSuing THE BRILLIANT QUALITY & CLIMBING THE PEAK OF THE SAME TRADE

ESWL ----LITHOTRIPTOR SET

800001

PNEUMATIC LITHOTRIPTOR



Speci:

Start & Stop pressure output
Pneumatic control parts made from Alloy material
Atraumatic

Functions:

Manual frequency adjustment
Manual air pressure adjustment

Technical Data:

Pressure output adjustment: 0--7 Bar
Pressure input: min 4 Bar
Pressure output: 1.5 -- 4 Bar (recommeend)
Power: 20Watts
Power supply: AC 100V -240 V 50/60 Hz
Stone smashing frequency: mono; 1 Hz ; 5 Hz ; 12 Hz

Dimension:

400mm x 112mm x 300mm (WxHxD)

Weight : 10 KGS

800002

COMPRESSED AIR SUPPLY



Speci:

Oil Free Self Lubrication air compressor
Two Cylinder
Low noise

Technical Data:

Power:750 Watts
Air flow rate: min 70l/ Minutes
Max. Pressure : 0.8 Mpa ±10%
Capacity of air bottle :38 L
Power Supply : AC 200V -240V
Noise: max 55dB

Dimention: 540mm x560mm x 1000mm (W x H x D)

Weight: 60 KGs

800003

ESWL HAND PIECE

Durable hand piece for ESWL design



PRODE

800004

φ 0.8 × 600

800005

φ 1.0 × 600

800006

φ 1.6 × 600

800007

φ 2.0 × 400

800008

SUCTION HAND PIECE



ESSE 3 ENDOSCOPE

PURSUING THE BRILLIANT QUALITY & CLIMBING THE PEAK OF THE SAME TRADE

ESWL ----LITHOTRIPTOR SET

800009

TUBE



800010

HANDLE



800011

HANDLE PIECE GUIDE



800012

FOOT SWITCH (WATERPROOF)



800013

SILICON TUBE



800014

PROTECTION CAP FOR STERILIZATION

