

RAIL GROUP

Model **GACB 2**
(code: 0500202)

It is hung on an universal rail, 30 x 10 mm. It is connected to the central vacuum system by a vacuum regulator to get a gradual suction and to use the system better.



- **Vacuum regulator mod. RV 1000 CE 0546** fitted with **ON-OFF** valve and vacuum gauge with scale -1000 to 0 mbar.
- **Chrome-plated universal bracket** for rail 30 x 10
- **Antibacterial filter** Ø 50 mm
- **Expansion chamber** for protecting the suction circuit, equipped with overfill valve (quick conical clutch)
- Connecting **clutch** to the central vacuum system
- **Bracket** and supporting **ring** for container
- **1 or 2 l. makralon collection container** which can be sterilized.
- **Silicone tube** Ø 7x12 mm
- **Single-patient collection tube**
- **Double funnel connection**
- **Dimensions:** cm. 30x11x40 h
- **Weight:** Kg. 3

OPERATING MODE: CONTINUOUS
POWER SOURCE: CENTRAL
VACUUM SYSTEM

IN CONFORMITY WITH:
UNI EN ISO 10079-3
93/42 EEC
VACUUM REGULATOR
RV 1000



0546 CLASS II a



ESSE 3 Via Garibaldi 30
14002 Castelnuovo D.B.(AT) – I
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

