

MEDICAL DENTAL AIR COMPRESSORS **ES**

ES 35 / HK-2EW-35

POWER SUPPLY : 220/240V , 50/60 HZ
CONSUMPTION : 3,7A (750W)
FLOW : 100 lt/MIN
TANK CAPACITY : 35 LT
OUTPUT : 3,5 BAR
NOISE : less of 60 dB
DICHARGE CAPACITY : 0,8 MPa



ES 45 / HK-3EW-45

POWER SUPPLY : 220/240V , 50/60 HZ
CONSUMPTION : 4,4A (1100W)
FLOW : 140 lt/MIN
TANK CAPACITY : 45 LT
OUTPUT : 3,5 BAR
NOISE : less of 60 dB
DICHARGE CAPACITY : 0,8 MPa



ESSE 3 Via Garibaldi 30
14002 Castelnuovo D.B.(AT) - I
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com



SILENT TOP AIR COMPRESSORS



SPECIFICATIONS

ES 30 - HK5.1.W30

Supply voltage: 220V AC 50Hz
 Electric current: 2.4A
 Frequency: 550W
 Power: 70L/min
 Capacity discharge: 0.8MPa
 Air container volume: 30L
 Noise: 56~65dB



CST192/60



CST192/100



CST193/130



CST193/200

CST193/240



Model	Power W	Discharge L/min	Max Pressure Mpa (4 bar)	Tank Capacity L	Noise dB(A)
CST192/60	550	60	0.8	25	52
CST192/100	1100	120	0.8	60	53
CST193/130	1500	160	0.8	60	54
CST193/200	1650	180	0.8	80	54
CST193/240	2250	240	0.8	80	55

Remark: It's available Noiseless box and Dryer as option.

We reserve the right to make any technical alterations without prior notice.