



OXYGEN GENERATOR

ΟΧΛΓΕΝ ΓΕΝΕΡΑΤΟΡ



EM-ZY2L



EM-ZY3LW



EM-ZY4LM



EM-ZY5LW

Specification	Code	EM-ZY2L	EM-ZY3LW	EM-ZY4LW	EM-ZY5LW
	Main voltage (V/Hz)		AC220/50	AC220/50	AC220/50
Power input (VA)		300	400	400	450
Oxygen content (%)		≥90	≥90	≥90	≥90
Rated oxygen generating efficiency (L/min)		2	3	4	5
Atomization efficiency (ml/min)		—	≥0.2	≥0.2	≥0.2
Mpa Oxygen delivery pressure (Mpa)		0.04-0.06	0.04-0.06	0.04-0.06	0.04-0.06
Noise level [dB(A)]		≤48	≤50	≤50	≤50
Relative humidity of operating environment (%)		≤80	≤80	≤80	≤80
Ambient temperature (°C)		5-40	5-40	5-40	5-40
Atmospheric Alarm for (Kpa)		86-106	86-106	86-106	86-106
Dimensions (mm)		270x430x475	300x560x575	300x560x575	300x560x575
Net weight (Kg)		20	23	25	30
Dual flow meter		—	—	—	✓
Alarm for power breatkoff		✓	✓	✓	✓
Alarm for excessive pressure		✓	✓	✓	✓
Alarm for adequate pressure		✓	✓	✓	✓
Atomizing device		—	—	✓	✓
Water drainage device		—	—	—	✓
Nebulization Function		—	✓	✓	✓
Fitering device		Round-shaped filter with tubular filtering element inside			Round-shaped primary filter Quadranglar secondary filter



Esse3 Via Garibaldi, 30 14022 Castelnuovo Don Bosco (AT) Italy

TEL: +39 011 992 77 06

FAX: +39 011 992 75 06

HTTP: WWW.4ci.com

E-MAIL: esse3@chierinet.it

