

ESSE 3

Trolley-9020

Endoscope use trolley



Dimensions:

LxWxH :
56cmx70cmx112cm (+13 Arm)

Draw height : 14 cm

Four Floor
(height from up to the bottom) :
220-200-200-240

Four ABS stable wheels
(with lock function)
Weight(Net): 45 KG



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

