



ES-DIAL

Electric dialyse chair



HEMODIALYSIS

ES-DIAL

... 9`YVWf]WX]U`ngY`W U]f



Cod. **ES-8=5@**

Comfortable:

- Wide adjustment range for the most comfortable position during dialysis treatment.
- Doctor or nurse can adjust the chair to the most optimal operation position according to the treatment requirement.
- The shock body-position can supply a perfect treatment condition for the patients with low blood pressure shock reactions.
- The real-leather soft mattress can slow down the body pressure point, it makes sure the blood vein of body will not be pressed at any situation.

Safe and Humanized design

- Anti-eroded, harmless and tire proof material, no sharp edge, very easy to be disinfected and cleaned.

More functions (optional)

- Lifter rank and shelf to supply the convenience of reading, writing, using computers or taking dinner
- USB port can supply more entertainment programs to play music and broadcast MP3 for patients
- Can be linked with hook and transfusion rack
- Weight scale function

Technical Features

Power supply	AC(220/110V) / (50/60Hz)	
Power	230 VA	
Insurance	1 A	
Noise	≤ 65 dB	
Weight	70 Kg	
Max. support weight	250 Kg	
DC 24V motor qty	3 pcs	
Safety type	II type	
Protection type	B type	
Colour	Rack	Light grey / White
	Mattress	Light blue / Light yellow / Light grey
Adjustable angle for back	Max. 75° ~ 80°	
Adjustable angle for small leg	Max. 30° ~ 40°	
Mattress thickness	5 ~ 7 cm	
Mattress size for hip	85 ~ 66 cm	
Mattress size for small leg	30 ~ 66 cm	
Seat cushion	58 ~ 66 cm	
Frame	Sponge with the real leather surface	
Complete size (LxWxH)	185 x 72 x 55 cm	



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

