



ESMT-1000

All-in-one Medical Computer Trolley



Specification

PC Specifications	Operate system	Windows 10 professional, 64 bit
	CPU	Intel® Celeron® CPU 3865U@1.8GHz
	RAM	4 GB
	Hard Drive	1 T
	Optical Drive	DVD-RW
	Input Devices	Alphanumeric Keyboard Wireless Mouse
	External Interface	5 USB
		1 RJ45 Network Ports 1 RS232 Serial Port
Power Supply	110V 2.0A / 220V 1.0A 50 Hz / 60 Hz	
Physical Specifications	Display	21.5" IPS-LCD (1920 × 1080)
	Dimensions	640 × 510 × 1235 mm (length × width × high)
	Weight	≤70 kg (with safe working load)
	Features	Wire accessory basket Printer tray Drawer for paper 4 trolley casters with lock pedal BP monitor fixing plate
	Load capacity	Basket: 5 kg Printer tray: 7.5 kg Drawer: 3.5 kg



Safety Specifications	Standards Compliance	IEC 60601-1:2005/A1:2012 EN 60601-1:2006/A1:2013 IEC 60601-1-2:2014 EN 60601-1-2:2015
	Anti-electric Shock Type	Class I equipment
	EMC	CISPR 11 Group 1 Class A
Environment Specifications	Temperature	Transport and Storage: -20 °C ~ +55 °C (-4°F ~ +131 °F) Working: +5 °C ~ + 40 °C (+41 °F ~ +104 °F)
	Relative Humidity	Transport and Storage: 25% RH ~ 93% RH (non-condensing) Working: 25% RH ~ 80% RH (non-condensing)
	Atmospheric Pressure	Transport and Storage: 70 kPa ~ 106 kPa Working: 86 kPa ~ 106 kPa



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

