



# ESE-FT100

Liver Diagnosis Ultrasound System



IMAGING SYSTEMS



*Shear Wave Quantificational Ultrasound Diagnostic System.  
An Aid to Diagnosis of Liver Fibrosis, Cirrhosis and Steatosis.*

*It facilitates early screening, early prevention, early diagnosis, and early treatment for chronic liver diseases.*

**Applied Technology:** *Transient Elastography*

**Light-weight:** *easy to carry and transport*

**Rapid:** *Fast Measurement & Immediate Results*

**Non invasive:** *Regular Disease Monitoring And Patient Follow-up*


## **Dynamic Wideband Fibrosis Scanning Probe**

- Applicable to different body sizes (adulte, overweight/obese patients, children)
- Highly succesfull detection rate
- Controlled probe pressure ensures high detection accuracy
- Optimized scientific algorithms provide fast and quantitative results



## **Innovative Technology for Liver Diagnosis**

- Non-Invasive • Rapid
- Reproducible • Accurate
- Quantitative • Reliable



### **Standard configuration:**

- Main unit with integrated battery
- 12.3" high-resolution touchscreen LCD monitor
- Fibrosis scanning probe
- Probe holder
- Foot switch

### **Option:**

- Wireless keyboard and mouse
- Customized trolley



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

