



# ESFD-806

## Fetal monitor




INTENSIVE CARES



# ESFD-806


## Fetal monitor



- 10.2" LED Touch Screen
- Screen rotation
- Mass storage device
- USB memory
- Waterproof probe
- Auto Interpretation



Touch Mode



Quick Input



Convenient Setting



Unlimited Storage

## Technical Specifications

- 10.2" LED touch-screen, 90° rotation screen.
- Twins Probe Available.
- CTG waveforms can be frozen.
- Thermal printer with 112mm, synchronization print.
- Auto interpretation and NST Grade.
- High performance waterproof Probe is equipped.
- Up to 500 patients data saving, USB can be supported.
- Data transmitting to PC via USB.
- Patient's Information can be Input.
- Historical data can be reviewed and printed.
- Built-in rechargeable li-ion battery, AC/DC power conversion.
- Sleep mode to save energy and extend LED Screen life.
- Portable design with hidden handle, convenient to carry.
- Multiple language support.
- 7 level printing thick can be chosen.
- Adapt to 110-230V, 50/60Hz AC power supply.

<b>Ultrasonic Operating Frequency</b>	3MHz, ±10%
<b>Systolic pressure</b>	0-100 unit
<b>Operation mode</b>	Manual/Auto
<b>Specification of fuse</b>	T1.6A 250V ø5x20
<b>Recording paper</b>	112x100mm Z-fold thermal paper
<b>FHR range</b>	50-240 beats/min
<b>Paper speed</b>	1, 2, 3 cm/min (±3%)
<b>Recorder</b>	High resolution thermal printer
<b>Battery</b>	Built-in rechargeable Li-ion battery 11.1V (2200mAh)
<b>Machine size</b>	307x255x91mm



ESSE3 srl, Via Garibaldi 30  
 14022 Castelnuovo D.B. (AT)  
 Tel +39 011 99 27 706  
 Fax +39 011 99 27 506  
 e-mail: [danielec@esse3-medical.com](mailto:danielec@esse3-medical.com)  
 web: [www.esse3-medical.com](http://www.esse3-medical.com)

