



ESUV-02S

UV-Transilluminator



LABORATORY



Features

UV Transilluminator is the necessary instrument for observed, photographed and cutting the electrophoresis gel.

- * With the operating port for cutting the gel.
- * No need for dark room, can be used for all weather.
- * Using the UV filter glass quartz ultraviolet lamp, long service life, light uniform.

Technical parameters

Transmission Size(W*L)	200*150mm
Transmission UV Wavelength	302nm
Reflectance UV Wavelength	254, 365nm
Transmission UV Lamp Power	302nm(8W)
Reflectance UV Lamp Power	254nm(11W), 365nm(11W)
Power Supply	AC110/220V±10%, 50/60Hz
External Size(W*D*H)	435*295*490 mm
Package Size(W*D*H)	700*640*810 mm
Net Weight	17kg
Gross Weight	45 kg



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

