



ES-DF 7

8 8 [] hU' < chd` UHY



LABORATORY

ES-DF)) 7

.....@98`8][]HJ`< chd`UH Y



Cod. ES-PR55C

- Large square plate
- Solid and sturdy working plate with a max temperature of 580 °C
- Temperature displayed on the LED digital display, easy to read
- Separated safety circuits with fixed safety temperature of 580°C
- Optional PT1000 temperature sensor
- "HOT" plate warning signal even if switched off to prevent accidental burns

Technical Features:

Dimension of working plate	184x184 mm
Working plate material	Ceramic
Power	1010 W
Voltage	100-120/200-240 VAC
Frequency	50/60
Heating power	1000 W
Heating Range	Up to 500°C, ±5 °C
Control accuracy of work plate	± 10 °C
Safety temperature	580°C
Temperature display	LED
Display temperature accuracy	± 1 °C
Display temperature accuracy with optional external temperature sensor PT1000	± 0.5 °C
Heating warning signal	50 °C
Protection class according to DIN EN60529	IP21
Dimensions	215x360x112 mm
Weight	4.5 Kg
Ambient temperature	5-40

Optional Accessories:

Cod. ES-C20100	Support clamp for PT-1000 sensor
Cod. ES-C20110	PT-1000 temperature sensor

Cod. **ES-PR55C** +
Cod. **ES-C20100** +
Cod. **ES-C20110**

