



**DATA SHEET** 

**Model: 318xx Line**

**Modular diagnostic station**






















The modular diagnostic station can be individually combined according to requirements and the room situation, it is variably extendable to the right and the left.

The basic module is available with one handle and can be supplemented with up to 4 further handle modules.

....

- 3.5V Xenon or LED illumination creates lighting conditions similar to daylight
- optimum electronic efficiency using an adjustable dimmer
- constant output voltage independent from the mains supply prolongs the service life of the lamps and guarantees optimum light intensity
- automatic switching on/off of handles when removing or replacing them in the handle unit saves energy and increases the service life of the lamps
- non-slip, ergonomically designed handles with rheostat: infinitely variable regulation of brightness as well as manual switching on/off possible on the handle
- spiral cord expandable up to 3 m
- developed and manufactured in Germany

Model	Description
<b>MODULAR DIAGNOSTIC STATION - XENON</b>	
<p><b>31800</b></p> 	<p>Diagnostic station 3.5-230V* 50/60 Hz-standard (includes 2 handles, Xenon oto and ophthalmo heads, ear specula dispenser and clock)</p>
<p><b>31805</b></p> 	<p>Diagnostic station 3.5-230V* 50/60 Hz-large (includes 2 handles, Xenon oto and ophthalmo heads, ear specula dispenser, clock and sphygmomanometer)</p>
<b>MODULAR DIAGNOSTIC STATION - LED</b>	
<p><b>31801</b></p> 	<p>Diagnostic station 3.5-230V* 50/60 Hz-standard (includes 2 handles, LED oto and ophthalmo heads, ear specula dispenser and clock)</p>
<p><b>31806</b></p> 	<p>Diagnostic station 3.5-230V* 50/60 Hz-large (includes includes 2 handles, LED oto and ophthalmo heads, ear specula dispenser, clock and sphygmomanometer)</p>
<b>RI-FORMER MODULES AND HEADS</b>	
<p><b>31809</b></p> 	<p>Basic diagnostic station 1 handle 3.5 - 230V* 50/60 Hz</p>
<p><b>31810</b></p> 	<p>Extension module 1 handle</p>
<p><b>31811</b></p> 	<p>L2 otoscope head - Led 3.5V</p>
<p><b>31812</b></p> 	<p>L1 otoscope head - Xenon 3.5V</p>
<p><b>31813</b></p> 	<p>L1 ophthalmoscope head - Xenon 3.5V</p>
<p><b>31814</b></p> 	<p>Retinoscope head - Xenon 3.5V</p>

<p><b>31815</b></p> 	Dermatoscope head - Led 3.5V
<p><b>31817</b></p> 	L1, L2 ear speculum Ø 2.5 mm** - autoclavable 134°C, Box/10
<p><b>31819</b></p> 	L1, L2 ear speculum Ø 4 mm** - autoclavable 134°C, Box/10
<p><b>31820</b></p> 	Ear specula dispenser - specula not included
<p><b>31822</b></p> 	Sphygmomanometer with basket
<p><b>31823</b></p> 	N infrared thermometer with 40 covers
<p><b>31824</b></p> 	Probe cover for 31823 - box of 100 - spare
<p><b>31825</b></p> 	Xenon L1 otoscope bulb 3.5V - spare
<p><b>31826</b></p> 	Xenon ophthalmoscope L1, L2, L3 bulb 3.5V - spare

\*available on request 120V - 50/60 Hz