



F-202

Ultrasonic Nebulizer
for Aerosol



OPERATING THEATRE



Description

- Excellent therapeutic performances coupled with an easy use.
Fast nebulization, low noise level, excellent micronisation of medication,
- Dimensions: 24,5 x 10 x 19 cm

Technical Specification

Weight Poids 600gr
Noise level \leq 55 dBA
Power supply 230V - 50/60 Hz
Medication cup capacity 12 ml
Nebulization up to 2,0 ml/min
Particles 0,5 - 5 micron
Frequency 1,7 MHz

Accessories standard

05.5012/1 Mask adult
05.5013/1 Mask child
05.5021/1 Mouth piece
05.5018.01 Corrugated tube
05.5030.01 Nebulizer bulb
ON REQUEST
05.5009 Nebulizer set complete with child face mask
05.5010 Nebulizer set complete with adult face mask



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

