



ES-SM/228

Plaster Saw



OPERATING THEATRE



First quality plaster saw. Ergonomic shape made of sturdy and non-slip plastic material, for synthetic and traditional plaster, available electric code ES-SM/228 or battery operated ES-SM/228E. The engine is electronically controlled to preserve its life and the gradual starting system also guarantees the mechanical components against excessive wear.

The oscillations of the blades can be adjusted to pleasure by a command placed on the handle.

ES-SM/228

Technical features of the ES-SM/228 electric plaster saw:

Power supply: 220-240 Volt 50 / 60Hz mains

Manual absorption: 250 Watt

Oscillations: 12,000 - 21,000 min-1

Weight: 1.2 kg without cable - 1.6 kg with cable

Cable length: 4.5 m

Noise: 68 - 82 dB (A)

Package contents: plaster saw, 1 blade for synthetic plaster diam. 50mm (code ES-SM/228G)
n.1 synthetic plaster blade diam. 65 mm code ES-SM/228S), n. 2 11mm keys.

ES-SM/228E

Technical characteristics of the ES-SM/228E plaster saw with rechargeable battery complete with battery charger:

Power supply: 220-240 Volt 50 / 60Hz mains battery charger

Oscillations: 11,000 - 18,000 min-1

Weight: 0.97 kg without Accu - 1.55 kg with Accu

Noise: 65 - 75 dB (A)

Charging time: approx. 40 min.

Package contents: plaster saw, 1 blade for synthetic plaster diam. 50mm (code ES-SM/228G)
n.1 synthetic plaster blade diam. 65 mm code ES-SM/228S), n. 2 11mm keys, charger and batteries.

ACCESSORIES

For ES-SM/228:

ES-SM/228P Traditional plaster blade 65 mm

ES-SM/228R Traditional plaster blade 50 mm

ES-SM/228S Synthetic plaster blade 65 mm

ES-SM/228G Synthetic plaster blade 50 mm



For ES-SM/228E

ES-SM/228P Traditional plaster blade 65 mm

ES-SM/228R Traditional plaster blade 50 mm

ES-SM/228S Synthetic plaster blade 65 mm

ES-SM/228G Synthetic plaster blade 50 mm

ES-SM/228T Accu Li-Ionen Battery 14.4 V / 4000 mAh

ES-SM/228U Battery charger: 220 - 240V / 110W



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

