



**DATA SHEET**



**Model: ESLD-200**

**Mobile Led Lamp**



### Description

1/ illumination	40000
2/ color temperature	4000 ± 500K
3/ light adjustable steps	----
4/ restore color index	95-98
5/ shadowless proportion	98%
6/ LED using life	>50000h
7/ up doctor head temperature	<1°C
8/ power	220V ± 22V 50Hz ± 1Hz
9/ package	carton
10/ packing size	117*51*29cm
11/ gross weight	45KG



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

