



304 MEDICAL GRADE SS

ES-MATRIX ECHO

(ELECTRO HYDRAULIC OT TABLE WITH MANUAL OVERRIDE SYSTEM)



SPECIAL FEATURES

Engineered for increasingly complex surgical scenes, **ES-MATRIX ECHO** comes with extreme loads in any position. The modular structure enables it to be adapted to match the size of the patient and in no time at all it can be equipped for any surgical discipline. With the characteristics of high automation low noise and high reliability, **ES-MATRIX ECHO** can be regarded as an ideal electro-hydraulic operating table.

- Modular design tabletop tailored for various surgical discipline needs. The easy click new mounting point enables the table adapted to wide variety of surgery and patient's height, waterproof design Mattress.
- Extreme weight load capacity of 300 kg with highest safety and stability.
- Standard override panel control makes the control safer.

TECHNICAL DATA

Length of Table	2050 mm
Width of Table Top	550 mm
Height (without Mattress)	Min. 750 to max. 1000 mm
Trendelenburg	+/- 25°
Lateral Tilt	+/- 20°
Back Section	+ /- 70°
Head Section	+/- 60°
Leg Section	0°-90°
Floor Locking System	Center Floor Locking
Top Sliding	5" Reverse & 5" Forward
Weight Bearing Capacity	300 kg

SALIENT FEATURES

- Split Leg Section
- Stainless Fitting
- Detachable, head, leg & pelvic section
- In built kidney bridge
- Five sectional radio – Translucent Top
- Base & Cover made of attractive, easy to clean Stainless Steel
- Sophisticated mechanics provide smooth, step less articulation of the table top for precise patient positioning
- Standby Battery Back-up with recharging circuit

TABLE POSITIONS

- Height Adjustment – Up & Down
- Lateral Tilt – Left & Right
- Trendelenburg & Reverse Trendelenburg
- Flex & Reflex Position
- Top Sliding



TABLE POSITIONS



Trendelenburg



Reverse Trendelenburg



Flex Position



Reflex Position



Lateral Tilt - Left & Right



Height Adjustment - Up & Down



Gynec Position



Kidney Position



STANDARD ACCESSORIES




ARM REST



SHOULDER SUPPORT



ANESTHETIC SCREEN



CLAMPS



LEG SUPPORT



MATTRESS



COMPATIBLE FOR

- Albee's Type Ortho Attachment
- Hanging Ortho Attachment
- Neuro Attachment

ES-KM1206/A

(ALBEE'S TYPE ORTHO ATTACHMENT)



- The Orthopedic Attachment can be used with any operating table
- In addition, a pelvic support plate is supported on top of the frame and is formed with an asymmetrical construction.
- Stirring Instruments allows flexion/raising of legs & stirring/string ball joints locked from foot end.
- Two traction devices on sliding spot with ball joint allowing flexions & rotations of foot with traction shoe with metallic sole
- Padded orthopedic knee crutches/provides support to the legs of the patient.

ES-KM1206/H

(HANGING TYPE ORTHO ATTACHMENT)

- The Orthopedic Attachment can be used with any operating table
- In addition, a pelvic support plate is supported on top of the frame and is formed with an asymmetrical construction.
- Stirring Instruments allows flexion/raising of legs & stirring/string ball joints locked from foot end.
- Two traction devices on sliding spot with ball joint allowing flexions & rotations of foot with traction shoe with metallic sole
- Padded orthopedic knee crutches/provides support to the legs of the patient.



STANDARD ACCESSORIES LIST

- Mattress
- Arm Rest
- Shoulder support
- Side support
- Anesthetic Screen
- Leg Support
- Clamps
- Remote
- Manual Override system