



# ES-ACP68

Auto Chart Projector



OPHTHALMOLOGY / ENT



## Specification

*There are: C landolt, E snellen, K LETTER, Children and Number.*

Projection distance	1.5m~ 6m
Projection magnification	30 x (at 5m)
Projection size	330mm(W) x270mm(H) (at 5m)
Chart	33 different pattern,60Pcs Masks.
Speed of chart conversion	One chart per 0.03s
Mask	1 open, 5 horizontal lines, 8 vertical lines, 21 single letters, 1 red/green
speed of mask conversion	One mask per 0.03s
Program	2 sets programs, each program contains up to 30 steps
Light source	LEDLamp
Auto-off function	After 10 minutes idle time.
Power source	AC 220V, 50Hz or 110V, 60Hz
Power consumption	80W
Dimensions	300mm(L) x 230mm(W) x 240mm(H)
Net weight	6Kg
Accessories	Remote control, Polarized metal screen, halogen lamp, polarized glasses, fuses(2), batteries(2)
Optional accessories	Floor stand, wall stand



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

